



model 7516BOX

eye wash bottle box

FEATURES & BENEFITS

OPTIONS

- ❑ 500 ml eye wash bottle: 500 ml eye wash bottle with 0.9 % sterile sodium chloride solution for rinsing away foreign bodies from the eye (dirt, dust etc.). Rinsing time approx. 5 min. Can also be used for refilling the eye wash stations. According to DIN EN 15154-4.
- ❑ 500 ml eye wash bottle: Eye wash station with 1 x 500 ml eye wash and pictogram. Rinsing time approx. 5 min.

For more information, visit www.hawesco.com or call (888) 640-4297.



SPECIFICATIONS

Emergency station in dust proof box with 1 x 200 ml pH Neutral + 1 x 500 ml eye wash. Especially suitable in places where there is a risk of eye injuries involving acids or alkali. Equipped with a mirror inside the lid. Rinsing time approx. 7 min.

APPLICATIONS

Emergency station in dust proof box with 1 x 200 ml pH Neutral + 1 x 500 ml eye wash. Especially suitable in places where there is a risk of eye injuries involving acids or alkali. Equipped with a mirror inside the lid. Rinsing time approx. 7 min.